

TABLE  $\Rightarrow$ 

Sr. No.	Volume of NaOH added	Conductance
	ml	$\Omega^{-1}$
1.	0	5.6
2.	0.5	5.0
3.	1.0	4.3
4.	1.5	3.5
5.	2.0	2.9
6.	2.5	2.4
7.	3.0	2.2
8.	3.5	2.4
9.	4.0	2.5
10.	4.5	2.7
11.	5.0	2.8
12.	5.5	3.0
13.	6.0	3.2
14.	6.5	3.6
15.	7.0	4.0
16.	7.5	4.5
17.	8.0	4.9
18.	8.5	5.4
19.	9.0	5.8
20.	9.5	6.3
21.	10.0	6.7
22.	10.5	7.1

## DETERMINATION OF THE STRENGTH OF A MIXTURE OF ACETIC ACID AND HYDROCHLORIC ACID BY CONDUCTOMETRY.

### AIM $\implies$

To estimate the strength of the mixture of acetic acid and hydrochloric acid present in a given mixture by conductometry.

### APPARATUS REQUIRED $\implies$

Conductivity bridge, conductivity cell, beaker, pipette, micro-burette, glass rod, etc.

### REAGENTS $\implies$

- 1) Hydrochloric acid (HCl) - 0.1N.
- 2) Acetic acid ( $\text{CH}_3\text{COOH}$ ) - 0.1N.
- 3) Sodium Hydroxide (NaOH) - 0.5N.

### PRINCIPLE $\implies$

The conductivity of the solution is related to the mobility of ions which in turn is related with the size of ions. When a mixture of acids like a strong acid (HCl) and weak acid ( $\text{CH}_3\text{COOH}$ ) is titrated against a strong base (NaOH). HCl reacts first followed by  $\text{CH}_3\text{COOH}$ . When titration of strong acid and strong base is carried out, there is a decrease in conductivity, as highly mobile  $\text{H}^+$  ions are replaced by  $\text{Na}^+$ .



23.	11.0	7.5
24.	11.5	8.0
25.	12.0	8.4
26.	12.5	8.8
27.	13.0	9.2
28.	13.5	9.6
29.	14.0	10.0
30.	14.5	10.4
31.	15.0	10.7

#### CALCULATIONS $\Rightarrow$

1) Strength of HCl.

Volume of mixture ( $V_1$ ) = 30 ml

Normality of NaOH ( $N_2$ ) = 0.5 N

Volume of NaOH ( $V_2$ ) = 2.7 ml

$$\therefore \text{Normality of HCl } (N_1) = \frac{V_2 \times N_2}{V_1} = \frac{2.7 \times 0.5}{20}$$

$$= \frac{1.35}{20} = 0.0675 \text{ N}$$

2) Strength of  $\text{CH}_3\text{COOH}$ .

Volume of mixture ( $V_1$ ) = 20 ml

Normality of NaOH ( $N_1$ ) = 0.5 N


Volume of NaOH ( $V_2 - V_1$ ) = 3.3 ml

$$\therefore \text{Normality of } \text{CH}_3\text{COOH } (N_2) = \frac{N_1 \times (V_2 - V_1)}{V_1}$$

$$= \frac{0.5 \times 3.3}{20} = \frac{1.65}{20} = 0.0825 \text{ N}$$



When the whole strong acid is consumed, base reacts with weak acid and conductivity increases as unionized weak acid because ionized salt.



After both the acids are consumed, there is a step increase in conductivity which gain endpoint and this increase in conductivity is due to fast moving  $\text{OH}^-$  ions from NaOH solution.

By the amount of base consumed, amount of acid present is calculated.

#### PROCEDURE $\implies$

- 1) The solution is diluted to 100ml, 20ml of which is pipetted out into a clean beaker and 100ml of distilled water is added.
- 2) Conductivity cell is dipped into the test solution and titrated against NaOH with strong conductance is measured for each addition of 0.5ml NaOH.
- 3) After neutralization, amount of acid present is determined by amount of NaOH consumption of strong acid and weak acid.
- 4) Volume of base consumed for strong acid and weak acid are determined by plotting a base added first, endpoint corresponds to strong acid into the other is for weak point.

RESULTS  $\Rightarrow$ 

- 1) The strength of  $\text{HCl}$  present in given solution =  $0.0675 \text{ N}$ .
- 2) The strength of  $\text{CH}_3\text{COOH}$  present in given whole of the solution =  $0.0825 \text{ N}$ .



y-axis

Scale:

On x-axis  
1cm = 1ml

On y-axis  
1cm = 10<sup>-1</sup>

