1) Standarization of HCL.

Sr. No.	Vol. of Na <sub>2</sub> CO <sub>3</sub> (V <sub>2</sub> )	Burette reading I F		Concordant value (V.)	Indicator
	m1			m1	
1.	20.0	0	19.8	19.8	methyl
2.	20.0	0	19.8		orange

2) Estimation of Na2CO3 and NaOH in a given mixeture.

Sr. No.	Volume of the unknown solution (V2)				
	m1		ml	ml	
1.	20	0	25	35	
2.	20	0	25	35	
	Concordant va	lue	A = 25	8 = 35	

## CALCULATIONS ->

Titration 1: STANDARDIZATION OF HCL.

Volume of HC2 = 19.8 m2 (V, m1-end point) Volume of Na<sub>2</sub>CO<sub>3</sub> (V<sub>2</sub>) = 20 m1

Page No.

DETERMINATION	OF	NQ2CO3	AND	NaOH
MIXTURE	BY	TITRA	TION.	a vitile

AIM -

To determine the amount of Na<sub>2</sub>CO<sub>3</sub> and NaOH in a mixture using HC1 acid.

APPARATUS REQUIRED =

Burette, pipette, standard measuring flack, beaker, funnel, wash bottle, etc.

REAGENTS ->

- 1) std. Hydrochloric acid.
- 2) sodium carbonate solution (0.1 N)
- 3) sodium hydroxide
- 4) methyl orange, phenolphthalein.

PRINCIPLE ==>

When a known volume of the mixture is titrated with HC1 in presence of phenolphthalein, the acid reacts with all the sodium hydroxide and with only

half of the carbonate.

when the titration is continued with methyl orange indicator, the remaining half of CO3 ions will be neutralized with HC1 at the end point.

A = all hydroxeide ions + carbonate ions

8 = half the carbonate ions after phenolphthalein endpoint

2A = all carbonate ions

Teacher's Signature \_\_\_\_

Normality of Na<sub>2</sub>CO<sub>3</sub> (N<sub>2</sub>) = 
$$0.05$$
 N

Normality of HC2 (N<sub>1</sub>) =  $\frac{V_2 \times N_2}{V_1}$  =  $\frac{20 \times 0.05}{19.8}$ 

=  $0.0505$  N

24 Titration 2: ESTIMATION OF AMOUNT OF Na2CO3 AND NaOH.

· Estimation of amount of Na2CO3.

Here, A = 25, B = 35, C = B - A = 35 - 25 = 10. Volume of HCl (V<sub>i</sub>) =  $2C = 2 \times 10 = 20 \text{ ml}$ 

(where, 'c' is titre value after phenolphthalein endpoint. i.e. C = B - A).

> Normality of HC2 (N<sub>1</sub>) = 0.0505 N (From titration 1) Volume of mixture (V<sub>2</sub>) = 20 ml

Normality of Na<sub>2</sub>CO<sub>3</sub> (N<sub>2</sub>) =  $\frac{V_1 \times N_1}{V_2}$  =  $\frac{20 \times 0.0505}{20}$ 

= 0.0505 N

:. Amount of Na<sub>2</sub>CO<sub>3</sub> present in whole of given solution = N<sub>2</sub> × Eq. wt of Na<sub>2</sub>CO<sub>3</sub> (53)

= 0.26765 gllit.

· Estimation of amount of NaOH.

Volume of HC1 (Vi) = 15m2 (A-C)

Normality of HC2 (Ni) = 0.0505 N (From titration 1)

Volume of mixeture (V2) = 20 m2

Normality of NaOH (N2) =  $\frac{V_1 \times N_1}{V_2}$  =  $\frac{15 \times 0.0505}{20}$ 

= 0.037875 N

Exp	ot. No Page No
	A-B = all hydroxide ions.
	PROGEDURE ===
	Titration 1: STANDARDIZATION OF HCL.
1>	
	## ## ## ## ## ## ## ## ## ## ## ## ##
2}	Pipette solution = std. Na2CO3 (0.1N)
3}	Indicator = methyl orange (2-3 drops)
4}	End point = yellow to orange.
	Titration 2: ESTIMATION OF THE MIXTURE (Na2CO3 + NaOH).
1}	Burette solution = std. Hydrochloric acid.
2}	Pipette solution = 20ml of the given mixture.
3}	Indicator = phenolphthalein (2-3 drops)
4>	End point = pink to colourless Continue the titration
	by adding methyl orange (2-3 drops) and
	the end point is the colour change from
	yellow to orange.

Teacher's Signature \_\_\_\_\_

... Amount of NaOH present in whole of the given solution = N2 × Eq. wt. of NaOH (40)

$$= 0.037875 \times 40 = 1.515$$
10

= 0.1515 glit.

EXF	t. No Page No
17	Amount of Na <sub>2</sub> CO <sub>3</sub> present in the given solution = 0.26765 gllit.
2>	Amount of NaOH present in the given solution = 0.1515 gllit.
	Teacher's Signature